User and maintenance manual and CE certification

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HOODS  
Series models C1 / C2 / C3 - serial number:

**prodotto da/produced by:**  
IRIDE International srl  
Via S.Giuseppe 5, Azzano Decimo - PN  
Italia

**emesso da/issued by**  
IRIDE International srl  
Via S.Giuseppe 5, Azzano Decimo - PN  
Italia

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### FUME HOOD GENERAL DESCRIPTION

The fume hood comprises a welding set, doors or dresser, vertical sliding tempered crystal sash with aluminium frame, suction equipment, control panel and plug sockets. The filtering equipment comprises two metal pre-filters with a honeycomb pattern and it also hosts 2 activated carbon filters to be periodically replaced.\* Activated carbon is a carbon which, through specific techniques, acquires a large absorptive capacity. Absorption is a specific physical phenomenon which allows an absorbing body to fix on its surface another substance (adsorbate). The external metal case is made of galvanised steel plate, whereas stainless steel is used internally. The frame of the tempered crystal sash is made of anodised aluminium. The fume hood must be connected to a flue to allow the emission of the fumes produced by the hot plate and the cylinder furnace. The hole must be 2,200 mm high from the ground and the hole diameter must be 120 mm wide. Iride International s.r.l. is not responsible for any connection carried out with cables not complying with the minimum thickness required by the "ELECTRICAL PANEL" diagram.

\* the metal pre-filters must be washed at least every 600 working hours; the cleaning must be performed using a vapour cleaning machine. Likewise, the activated carbon filters must be replaced every 600 working hours, considering the average use of the hood by 6 users.

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## WIRING DIAGRAM COMPONENTS

- 1- residual current operated circuit-breaker (without over-current release) 2P 25A 0.03A "AC"
- 2- thermomagnetic switch 1P+N 16A "C"  
800 hood - 1 plug socket  
1,200 hood - 2 plug sockets  
1,600 hood - 3 plug sockets
- 3- 1P 10A 3x38 delayed fuse disconnecter door.
- 4- 2P 25A switch
- 5- digital time switch - weekly model
- 6- 2P 25A switch

## OTHER HOOD COMPONENTS

- switch 1, mobile CE plug 2P+T 16A 220V  
switch 2, wall CE plug sockets 2P+T 16A 220V  
switch 3, 60W oval overhead light with CAGE

## SYSTEM DESCRIPTION

The electrical system changes according to the model required (800, 1,200 or 1,600 hood).

The 800 hood uses cables with a carrying capacity of 4 mm<sup>2</sup>, therefore the maximum carrying capacity is 4,500W.


The 1,200 and 1,600 hoods use cables with a carrying capacity of 6 mm<sup>2</sup>, therefore the maximum carrying capacity is 6,000W. The maximum capacity is the electrical input of the entire hood. For example, the following hoods:

(800) 1 lamp 60W+motor 500W+plug 1-2,500W = 3,060W

(1200) 1 lamp 60W+motor 500W+plug 1-2,500W+plug 2-2,500W = 5,560W

(1600) 1 lamp 60W+motor 500W+plug 1-2,500W+plug 2-2,000W+plug 3-1,000W = 5,560W

**WARNING!** The total sum of the plugs must not exceed 6,000W, hence the various plugs may have different values but the above mentioned maximum value must never be exceeded.

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

*Within the watertight case there are 2 free clamps, brown and blue; in order to connect the ground, connect the copper rod. A spare fuse and clock instructions are also supplied with the package.*

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*This connection may be used with a Ø 80mm pipe communicating with the outside to allow the hermetically sealed hood to work with an internal air recycle.*

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## FILTERING

*This sophisticated heavy dust filtering system allows dusts to be conveyed into the hollow space on the right side of the hood.*

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
*To replace it, see figures no. 1, 2 and 3.*

*To reassemble, repeat the steps in reverse order.*

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## HERMETICAL SEALING

*Usually there are two clamps to prevent the glass from closing completely. When the user is following the instructions of item (A)-AIR INLET (recommended), the clamps can be removed to allow for the hermetical sealing of the hood and consequently to use the internal recycle system.*

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## WELDING SET

*Welding set, adjustable in height and a 360° revolving table.*

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**IRIDE INTERNATIONAL S.R.L.**

## DECLARATION OF CONFORMITY

We, IRIDE INTERNATIONAL S.r.l. – furniture for dentist offices and dental laboratories, with registered headquarters in Via San Giuseppe 5 - 33082 Azzano Decimo - Località Tiezzo - PN:

**HEREBY DECLARE UNDER OUR OWN RESPONSIBILITY THAT THE PRODUCT "HOODS" TO WHICH THIS DECLARATION RELATES WERE DESIGNED AND BUILT IN EXCLUSIVE MANUFACTURING IN CONFORMITY WITH THE FOLLOWING HARMONISED STANDARDS:**

**IEC EN 60601-1(IEC 65-5) of 1991 and its 1994 amendments**

**IEC EN 60204-1(IEC 44/5 and 44/6)**

**AND THE REQUIREMENTS OF THE DIRECTIVE:**

**2006/42/CE**

**2006/95/CE**

**2004/108/CE**

**IRIDE INTERNATIONAL S.R.L.**

**IRIDE INTERNATIONAL s.r.l.**

Via S. Giuseppe 5 | 33082 Azzano Decimo - PN

tel. +39 0434 633242 | fax +39 0434 632546

welcome@irideinternational.com | www.irideinternational.com